



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human FGF2/bFGF Antibody (SAA0437), PerCP

#### OVERVIEW

---

Reference	ARO-A10001
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

#### PRODUCT INFORMATION

---

Host Species	Mouse
Target	bFGF, FGF2, Basic fibroblast growth factor, HBGF-2, FGFB, FGF-2, Heparin-binding growth factor 2, Fibroblast growth factor 2
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-bFGF, FGF2, Basic fibroblast growth factor, HBGF-2, FGFB, FGF-2, Heparin-binding growth factor 2, Fibroblast growth factor 2