

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human FGF2/bFGF Antibody (SAA0437), PerCP
PRODUCT REFERENCE	ARO-A10001
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	bFGF, FGF2, Basic fibroblast growth factor, HBGF-2, FGFB, FGF-2, Heparin-binding growth factor 2, Fibroblast growth factor 2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 0.2% BSA, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-bFGF, FGF2, Basic fibroblast growth factor, HBGF-2, FGFB, FGF-2, Heparin-binding growth factor 2, Fibroblast growth factor 2