



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD41/ITGA2B & CD61/ITGB3 Antibody (OPG2), PerCP

#### OVERVIEW

|           |                     |
|-----------|---------------------|
| Reference | ARO-A10025          |
| Clonality | Monoclonal Antibody |
| Isotype   | IgG1, kappa         |

#### PRODUCT INFORMATION

|                     |                                                                                                                                                                                                                                                                                         |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host Species        | Mouse                                                                                                                                                                                                                                                                                   |
| Target              | CD41, ITGAB, GPIIb, ITGA2B, Platelet membrane glycoprotein IIb, GP2B, Integrin alpha-IIb, GPalpha IIb, GPIIIa, Integrin beta-3, CD61, ITGB3, Platelet membrane glycoprotein IIIa, GP3A                                                                                                  |
| Form                | Liquid                                                                                                                                                                                                                                                                                  |
| Buffer              | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.                                                                                                                                                                                                                                         |
| Stability & Storage | 4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| Aliases /Synonyms   | Anti-CD41, ITGAB, GPIIb, ITGA2B, Platelet membrane glycoprotein IIb, GP2B, Integrin alpha-IIb, GPalpha IIb, GPIIIa, Integrin beta-3, CD61, ITGB3, Platelet membrane glycoprotein IIIa, GP3A                                                                                             |