



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human PAC-1 Antibody (PAC-1), PerCP

OVERVIEW

| | |
|-----------|---------------------|
| Reference | ARO-A10026 |
| Clonality | Monoclonal Antibody |
| Isotype | IgM, kappa |

PRODUCT INFORMATION

| | |
|---------------------|---|
| Host Species | Mouse |
| Target | CD41, ITGAB, GPIIb, ITGA2B, Platelet membrane glycoprotein IIb, GP2B, Integrin alpha-IIb, GPalpha IIb, GPIIIa, Integrin beta-3, CD61, ITGB3, Platelet membrane glycoprotein IIIa, GP3A, PAC-1 |
| Form | Liquid |
| Buffer | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300. |
| Stability & Storage | 4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| Aliases /Synonyms | Anti-CD41, ITGAB, GPIIb, ITGA2B, Platelet membrane glycoprotein IIb, GP2B, Integrin alpha-IIb, GPalpha IIb, GPIIIa, Integrin beta-3, CD61, ITGB3, Platelet membrane glycoprotein IIIa, GP3A, PAC-1 |