



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD271/NGFR Antibody (T-2), PerCP

#### OVERVIEW

Reference	ARO-A10038
Clonality	Monoclonal Antibody
Isotype	IgG2a

#### PRODUCT INFORMATION

Host Species	Mouse
Target	Low affinity neurotrophin receptor p75NTR, TNFRSF16, CD271, Gp80-LNGFR, p75 ICD, Low-affinity nerve growth factor receptor, NGFR, NGF receptor, Tumor necrosis factor receptor superfamily member 16
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Low affinity neurotrophin receptor p75NTR, TNFRSF16, CD271, Gp80-LNGFR, p75 ICD, Low-affinity nerve growth factor receptor, NGFR, NGF receptor, Tumor necrosis factor receptor superfamily member 16