



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human PSA/KLK3 Antibody (SAA0506), PerCP

#### OVERVIEW

Reference	ARO-A10058
Clonality	Monoclonal Antibody
Isotype	IgG1

#### PRODUCT INFORMATION

Host Species	Human
Target	Prostate-specific antigen, KLK3, Kallikrein-3, APS, Semenogelase, PSA, Gamma-seminoprotein, Semin, P-30 antigen
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Prostate-specific antigen, KLK3, Kallikrein-3, APS, Semenogelase, PSA, Gamma-seminoprotein, Semin, P-30 antigen