



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human ALPP/PLAP Antibody (17E2#), PerCP

#### OVERVIEW

Reference	ARO-A10086
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	PLAP-1, ALPP, Placental alkaline phosphatase 1, PLAP, Alkaline phosphatase Regan isozyme, Alkaline phosphatase, placental type
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-PLAP-1, ALPP, Placental alkaline phosphatase 1, PLAP, Alkaline phosphatase Regan isozyme, Alkaline phosphatase, placental type