



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD18/ITGB2 Antibody (1B4), PerCP

OVERVIEW

Reference	ARO-A10091
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

PRODUCT INFORMATION

Host Species	Mouse
Target	Complement receptor C3 subunit beta, Cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta, Integrin beta-2, MFI7, CD18, ITGB2
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Complement receptor C3 subunit beta, Cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta, Integrin beta-2, MFI7, CD18, ITGB2