



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human IGF1 Antibody (SAA0432), PerCP

#### OVERVIEW

Reference	ARO-A10100
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	IGF-I, IGF1, Mechano growth factor, MGF, Insulin-like growth factor I, IBP1, Somatomedin-C
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-IGF-I, IGF1, Mechano growth factor, MGF, Insulin-like growth factor I, IBP1, Somatomedin-C