



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Mouse Thy-1.2/CD90.2 Antibody (30-H12), PerCP

#### OVERVIEW

Reference	ARO-A10112
Clonality	Monoclonal Antibody
Isotype	IgG2b, kappa

#### PRODUCT INFORMATION

Host Species	Rat
Target	THY1, CDw90, Thy-1 antigen, Thy-1 membrane glycoprotein, CD90
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-THY1, CDw90, Thy-1 antigen, Thy-1 membrane glycoprotein, CD90