



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD71/TFRC Antibody (SAA1396), PerCP

OVERVIEW

Reference	ARO-A10126
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

PRODUCT INFORMATION

Host Species	Mouse
Target	Transferrin receptor protein 1, TfR1, T9, Trfr, p90, CD71, TFRC, TfR, sTfR, TR
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Transferrin receptor protein 1, TfR1, T9, Trfr, p90, CD71, TFRC, TfR, sTfR, TR