



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CXCL10/IP-10 Antibody (SAA0413), PerCP

#### OVERVIEW

Reference	ARO-A10129
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

#### PRODUCT INFORMATION

Host Species	Rat
Target	Gamma-IP10, C-X-C motif chemokine 10, INP10, 10 kDa interferon gamma-induced protein, Small-inducible cytokine B10, SCYB10, IP-10, CXCL10
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Gamma-IP10, C-X-C motif chemokine 10, INP10, 10 kDa interferon gamma-induced protein, Small-inducible cytokine B10, SCYB10, IP-10, CXCL10