



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human PPBP/CXCL7/NAP-2 Antibody (SAA0464), PerCP

#### OVERVIEW

Reference	ARO-A10131
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	NAP-2, NAP-2(1-63), NAP-2(73), PPBP, CXCL7, PBP, Small-inducible cytokine B7, Macrophage-derived growth factor, LA-PF4, Leukocyte-derived growth factor, CTAP-III, CTAP-III(1-81), Low-affinity platelet factor IV, TGB1, Beta-TG, C-X-C motif chemokine 7, MDG
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-NAP-2, NAP-2(1-63), NAP-2(73), PPBP, CXCL7, PBP, Small-inducible cytokine B7, Macrophage-derived growth factor, LA-PF4, Leukocyte-derived growth factor, CTAP-III, CTAP-III(1-81), Low-affinity platelet factor IV, TGB1, Beta-TG, C-X-C motif chemokine 7