



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Mouse CD25/IL2RA Antibody (PC61/PC61.5.3), PerCP

#### OVERVIEW

Reference	ARO-A10166
Clonality	Monoclonal Antibody
Isotype	IgG1, lambda

#### PRODUCT INFORMATION

Host Species	Rat
Target	Interleukin-2 receptor subunit alpha, IL-2 receptor subunit alpha, IL-2-RA, IL-2R subunit alpha, IL2-RA, p55, CD25, IL2ra, IL2r
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Interleukin-2 receptor subunit alpha, IL-2 receptor subunit alpha, IL-2-RA, IL-2R subunit alpha, IL2-RA, p55, CD25, IL2ra, IL2r