



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human IL33 Antibody (SAA0395), PerCP

OVERVIEW

Reference	ARO-A10191
Clonality	Monoclonal Antibody
Isotype	IgG1, lambda

PRODUCT INFORMATION

Host Species	Mouse
Target	NF-HEV, IL-1F11, C9orf26, IL-33, Interleukin-1 family member 11, IL33, Nuclear factor from high endothelial venules, NFHEV, Interleukin-33, IL1F11
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-NF-HEV, IL-1F11, C9orf26, IL-33, Interleukin-1 family member 11, IL33, Nuclear factor from high endothelial venules, NFHEV, Interleukin-33, IL1F11