

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Anti-Human GDF11/BMP11 Antibody (SAA0431), PerCP |
| PRODUCT REFERENCE | ARO-A10197 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG2a, kappa |
| TARGET | GDF11, BMP-11, BMP11, Growth/differentiation factor 11, GDF-11, Bone morphogenetic protein 11 |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| ORIGIN SPECIES | <i>Mus musculus (Mouse)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 0.2% BSA, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-GDF11, BMP-11, BMP11, Growth/differentiation factor 11, GDF-11, Bone morphogenetic protein 11 |