



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human GDF11/BMP11 Antibody (SAA0431), PerCP

#### OVERVIEW

---

Reference	ARO-A10197
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

#### PRODUCT INFORMATION

---

Host Species	Mouse
Target	GDF11, BMP-11, BMP11, Growth/differentiation factor 11, GDF-11, Bone morphogenetic protein 11
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-GDF11, BMP-11, BMP11, Growth/differentiation factor 11, GDF-11, Bone morphogenetic protein 11