



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD218b/IL18RAP Antibody (3131#), PerCP

OVERVIEW

Reference	ARO-A10200
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

PRODUCT INFORMATION

Host Species	Mouse
Target	AcPL, Accessory protein-like, Interleukin-1 receptor 7, IL-18R-beta, CD218 antigen-like family member B, IL-1RAcPL, Interleukin-18 receptor accessory protein, CDw218b, IL-1R-7, IL1R7, IL-1R7, IL-18RAcP, CD218b, IL-1R accessory protein-like, IL18RAP, Inter
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-AcPL, Accessory protein-like, Interleukin-1 receptor 7, IL-18R-beta, CD218 antigen-like family member B, IL-1RAcPL, Interleukin-18 receptor accessory protein, CDw218b, IL-1R-7, IL1R7, IL-1R7, IL-18RAcP, CD218b, IL-1R accessory protein-like, IL18RAP,