



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD256/TNFSF13 Antibody (VIS-649), PerCP

#### OVERVIEW

Reference	ARO-A10206
Clonality	Monoclonal Antibody
Isotype	IgG2, Kappa

#### PRODUCT INFORMATION

Host Species	Human
Target	TNF- and APOL-related leukocyte expressed ligand 2, APRIL, CD256, Tumor necrosis factor ligand superfamily member 13, TNF-related death ligand 1, TALL-2, TNFSF13, ZTNF2, A proliferation-inducing ligand, TALL2, TRDL-1
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-TNF- and APOL-related leukocyte expressed ligand 2, APRIL, CD256, Tumor necrosis factor ligand superfamily member 13, TNF-related death ligand 1, TALL-2, TNFSF13, ZTNF2, A proliferation-inducing ligand, TALL2, TRDL-1