



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD358/TNFRSF21 Antibody (M53-E04), PerCP

OVERVIEW

Reference	ARO-A10207
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

PRODUCT INFORMATION

Host Species	Mouse
Target	DR6, Death receptor 6, Tumor necrosis factor receptor superfamily member 21, TNFRSF21, CD358
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-DR6, Death receptor 6, Tumor necrosis factor receptor superfamily member 21, TNFRSF21, CD358