



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD275/ICOSLG Antibody (SAA0089), PerCP

#### OVERVIEW

Reference	ARO-A10209
Clonality	Monoclonal Antibody
Isotype	IgG2, kappa

#### PRODUCT INFORMATION

Host Species	Human
Target	B7RP1, B7-related protein 1, CD275, ICOSL, ICOSLG, KIAA0653, B7 homolog 2, ICOS ligand, B7RP-1, B7H2, B7-H2, B7-like protein GI50
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-B7RP1, B7-related protein 1, CD275, ICOSL, ICOSLG, KIAA0653, B7 homolog 2, ICOS ligand, B7RP-1, B7H2, B7-H2, B7-like protein GI50