



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD16b/FCGR3B Antibody (B88-9), PerCP

OVERVIEW

Reference	ARO-A10211
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

PRODUCT INFORMATION

Host Species	Mouse
Target	IGFR3, FCGR3B, CD16b, Low affinity immunoglobulin gamma Fc region receptor III-B, Fc-gamma RIII, IgG Fc receptor III-1, FCGR3, CD16B, Fc-gamma RIII-beta, Fc-gamma RIIIb, FCG3, FcRIII, FcR-10, FcRIIIb
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-IGFR3, FCGR3B, CD16b, Low affinity immunoglobulin gamma Fc region receptor III-B, Fc-gamma RIII, IgG Fc receptor III-1, FCGR3, CD16B, Fc-gamma RIII-beta, Fc-gamma RIIIb, FCG3, FcRIII, FcR-10, FcRIIIb