



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CXCL13/BCA-1/BLC Antibody (SAA0468), PerCP

OVERVIEW

Reference	ARO-A10221
Clonality	Monoclonal Antibody
Isotype	IgG2a, lambda

PRODUCT INFORMATION

Host Species	Mouse
Target	SCYB13, BLC, CXCL13, B cell-attracting chemokine 1, B lymphocyte chemoattractant, Small-inducible cytokine B13, CXC chemokine BLC, Angie, C-X-C motif chemokine 13, BCA-1, BCA1
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-SCYB13, BLC, CXCL13, B cell-attracting chemokine 1, B lymphocyte chemoattractant, Small-inducible cytokine B13, CXC chemokine BLC, Angie, C-X-C motif chemokine 13, BCA-1, BCA1