



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CLDN3/4 Antibody (KM3907), PerCP

#### OVERVIEW

Reference	ARO-A10226
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	Claudin-3, Rat ventral prostate.1 protein homolog, CPE-R 2, C7orf1, Clostridium perfringens enterotoxin receptor 2, CPE-receptor 2, CPETR2, hRVP1, CLDN3, CPER, CLDN4, Williams-Beuren syndrome chromosomal region 8 protein, CPE-R, CPE-receptor, WBSCR8, Clau
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Claudin-3, Rat ventral prostate.1 protein homolog, CPE-R 2, C7orf1, Clostridium perfringens enterotoxin receptor 2, CPE-receptor 2, CPETR2, hRVP1, CLDN3, CPER, CLDN4, Williams-Beuren syndrome chromosomal region 8 protein, CPE-R, CPE-receptor, WBSCR8,