



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD337/NCR3 Antibody (BGA1831), PerCP

OVERVIEW

Reference	ARO-A10229
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

PRODUCT INFORMATION

Host Species	Mouse
Target	NKp30, LY117, NCR3, 1C7, NK-p30, CD337, Natural cytotoxicity triggering receptor 3, Activating natural killer receptor p30, Natural killer cell p30-related protein
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-NKp30, LY117, NCR3, 1C7, NK-p30, CD337, Natural cytotoxicity triggering receptor 3, Activating natural killer receptor p30, Natural killer cell p30-related protein