



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD254/RANKL/TNFSF11 Antibody (SAA0080), PerCP

OVERVIEW

Reference	ARO-A10230
Clonality	Monoclonal Antibody
Isotype	IgG2, kappa

PRODUCT INFORMATION

Host Species	Human
Target	ODF, TNF-related activation-induced cytokine, Osteoclast differentiation factor, TNFSF11, TRANCE, Osteoprotegerin ligand, RANKL, Tumor necrosis factor ligand superfamily member 11, Receptor activator of nuclear factor kappa-B ligand, CD254, OPGL
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-ODF, TNF-related activation-induced cytokine, Osteoclast differentiation factor, TNFSF11, TRANCE, Osteoprotegerin ligand, RANKL, Tumor necrosis factor ligand superfamily member 11, Receptor activator of nuclear factor kappa-B ligand, CD254, OPGL