



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Mouse CD304/NRP1/NRP2 Antibody (YW68.11.26), PerCP

#### OVERVIEW

Reference	ARO-A10232
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

#### PRODUCT INFORMATION

Host Species	Rat
Target	NRP1, Neuropilin-1, VEGF165R, NRP, CD304, Vascular endothelial cell growth factor 165 receptor, Vascular endothelial cell growth factor 165 receptor 2, Neuropilin-2, VEGF165R2, NRP2
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-NRP1, Neuropilin-1, VEGF165R, NRP, CD304, Vascular endothelial cell growth factor 165 receptor, Vascular endothelial cell growth factor 165 receptor 2, Neuropilin-2, VEGF165R2, NRP2