



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD156c/ADAM10 Antibody (11G2), PerCP

#### OVERVIEW

Reference	ARO-A10239
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	ADAM 10, Mammalian disintegrin-metalloprotease, CD156c, MADM, ADAM10, Kuzbanian protein homolog, KUZ, CDw156, Disintegrin and metalloproteinase domain-containing protein 10
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-ADAM 10, Mammalian disintegrin-metalloprotease, CD156c, MADM, ADAM10, Kuzbanian protein homolog, KUZ, CDw156, Disintegrin and metalloproteinase domain-containing protein 10