



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD261/TNFRSF10A/DR4 Antibody (SAA0083), PerCP

#### OVERVIEW

|           |                     |
|-----------|---------------------|
| Reference | ARO-A10248          |
| Clonality | Monoclonal Antibody |
| Isotype   | IgG1                |

#### PRODUCT INFORMATION

|                     |   |
|---------------------|---|
| Host Species        | Human   |
| Target              | TNF-related apoptosis-inducing ligand receptor 1, APO2, Death receptor 4, TRAIL receptor 1, TRAIL-R1, CD261, DR4, TRAILR1, TNFRSF10A, Tumor necrosis factor receptor superfamily member 10A   |
| Form                | Liquid  |
| Buffer              | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.   |
| Stability & Storage | 4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| Aliases /Synonyms   | Anti-TNF-related apoptosis-inducing ligand receptor 1, APO2, Death receptor 4, TRAIL receptor 1, TRAIL-R1, CD261, DR4, TRAILR1, TNFRSF10A, Tumor necrosis factor receptor superfamily member 10A  |