



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CCL24/Eotaxin-2/MPIF-2 Antibody (SAA0461), PerCP

OVERVIEW

Reference	ARO-A10250
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

PRODUCT INFORMATION

Host Species	Mouse
Target	MPIF2, Eosinophil chemotactic protein 2, Small-inducible cytokine A24, CCL24, Eotaxin-2, SCYA24, MPIF-2, Myeloid progenitor inhibitory factor 2, CK-beta-6, C-C motif chemokine 24
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-MPIF2, Eosinophil chemotactic protein 2, Small-inducible cytokine A24, CCL24, Eotaxin-2, SCYA24, MPIF-2, Myeloid progenitor inhibitory factor 2, CK-beta-6, C-C motif chemokine 24