



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD126/IL6R/IL-6RA Antibody (SAA0048), APC

#### OVERVIEW

Reference	ARO-A10256
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	sIL6R, Interleukin-6 receptor subunit alpha, CD126, IL-6R 1, IL-6R subunit alpha, Membrane glycoprotein 80, gp80, IL6R, IL-6RA, IL-6 receptor subunit alpha, IL-6R-alpha
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-sIL6R, Interleukin-6 receptor subunit alpha, CD126, IL-6R 1, IL-6R subunit alpha, Membrane glycoprotein 80, gp80, IL6R, IL-6RA, IL-6 receptor subunit alpha, IL-6R-alpha