



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD23 Antibody (MHM6), APC

#### OVERVIEW

Reference	ARO-A10318
Clonality	Monoclonal Antibody
Isotype	IgG2a

#### PRODUCT INFORMATION

Host Species	Mouse
Target	Lymphocyte IgE receptor, C-type lectin domain family 4 member J, Low affinity immunoglobulin epsilon Fc receptor, CD23A, BLAST-2, Immunoglobulin E-binding factor, IGEBF, CD23, Fc-epsilon-RII, CLEC4J, FCE2, FCER2
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Lymphocyte IgE receptor, C-type lectin domain family 4 member J, Low affinity immunoglobulin epsilon Fc receptor, CD23A, BLAST-2, Immunoglobulin E-binding factor, IGEBF, CD23, Fc-epsilon-RII, CLEC4J, FCE2, FCER2