



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD18/ITGB2 & CD11/ITGAL Antibody (23F2G), APC

#### OVERVIEW

Reference	ARO-A10343
Clonality	Monoclonal Antibody
Isotype	IgG4

#### PRODUCT INFORMATION

Host Species	Human
Target	Complement receptor C3 subunit beta, Cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta, Integrin beta-2, MFI7, CD18, ITGB2, CD11A, Leukocyte function-associated molecule 1 alpha chain, Leukocyte adhesion glycoprotein LFA-1 alpha chain, LF
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Complement receptor C3 subunit beta, Cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta, Integrin beta-2, MFI7, CD18, ITGB2, CD11A, Leukocyte function-associated molecule 1 alpha chain, Leukocyte adhesion glycoprotein LFA-1 alpha chai