



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD230/PRNP Antibody (ICSM 18), APC

#### OVERVIEW

Reference	ARO-A10369
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	PRNP, PrP27-30, PRIP, PrP, CD230, ALTPRP, ASCR, PrP33-35C, PRP, Major prion protein
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-PRNP, PrP27-30, PRIP, PrP, CD230, ALTPRP, ASCR, PrP33-35C, PRP, Major prion protein