



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human IL1B/IL1F2 Antibody (SAA0369), APC

#### OVERVIEW

Reference	ARO-A10422
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	Interleukin-1 beta, IL1B, Catabolin, IL1F2, IL-1 beta,IL-1 $\hat{1}^2$
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Interleukin-1 beta, IL1B, Catabolin, IL1F2, IL-1 beta,IL-1 $\hat{1}^2$