



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human TGFA/TGF-alpha Antibody (SAA0442), APC

#### OVERVIEW

Reference	ARO-A10434
Clonality	Monoclonal Antibody
Isotype	IgG2a, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	ETGF, Protransforming growth factor alpha, TGF type 1, TGFA, EGF-like TGF, TGF-alpha
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-ETGF, Protransforming growth factor alpha, TGF type 1, TGFA, EGF-like TGF, TGF-alpha