

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD160 Antibody (H7), APC
PRODUCT REFERENCE	ARO-A10439
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Natural killer cell receptor BY55, CD160 antigen, CD160, BY55

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 0.2% BSA, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Natural killer cell receptor BY55, CD160 antigen, CD160, BY55