



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD85a/LILRB3 Antibody (A1), APC

#### OVERVIEW

Reference	ARO-A10461
Clonality	Monoclonal Antibody
Isotype	IgG1, lambda

#### PRODUCT INFORMATION

Host Species	Human
Target	CD85a, ILT-5, Monocyte inhibitory receptor HL9, LILRB3, Immunoglobulin-like transcript 5, CD85 antigen-like family member A, ILT5, Leukocyte immunoglobulin-like receptor 3, LIR-3, Leukocyte immunoglobulin-like receptor subfamily B member 3, LIR3
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-CD85a, ILT-5, Monocyte inhibitory receptor HL9, LILRB3, Immunoglobulin-like transcript 5, CD85 antigen-like family member A, ILT5, Leukocyte immunoglobulin-like receptor 3, LIR-3, Leukocyte immunoglobulin-like receptor subfamily B member 3, LIR3