



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD203c/ENPP3 Antibody (H1.8), APC

OVERVIEW

| | |
|-----------|---------------------|
| Reference | ARO-A10489 |
| Clonality | Monoclonal Antibody |
| Isotype | IgG1, kappa |

PRODUCT INFORMATION

| | |
|---------------------|---|
| Host Species | Mouse |
| Target | NPPase, PD-Ibeta, Phosphodiesterase I beta, PDNP3, CD203c, Ectonucleotide pyrophosphatase/phosphodiesterase family member 3, E-NPP 3, ENPP3, NPP3, Phosphodiesterase I/nucleotide pyrophosphatase 3, Nucleotide diphosphatase |
| Form | Liquid |
| Buffer | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300. |
| Stability & Storage | 4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| Aliases /Synonyms | Anti-NPPase, PD-Ibeta, Phosphodiesterase I beta, PDNP3, CD203c, Ectonucleotide pyrophosphatase/phosphodiesterase family member 3, E-NPP 3, ENPP3, NPP3, Phosphodiesterase I/nucleotide pyrophosphatase 3, Nucleotide diphosphatase |