



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD261/TNFRSF10A/DR4 Antibody (SAA0083), APC

OVERVIEW

Reference	ARO-A10497
Clonality	Monoclonal Antibody
Isotype	IgG1

PRODUCT INFORMATION

Host Species	Human
Target	TNF-related apoptosis-inducing ligand receptor 1, APO2, Death receptor 4, TRAIL receptor 1, TRAIL-R1, CD261, DR4, TRAILR1, TNFRSF10A, Tumor necrosis factor receptor superfamily member 10A
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-TNF-related apoptosis-inducing ligand receptor 1, APO2, Death receptor 4, TRAIL receptor 1, TRAIL-R1, CD261, DR4, TRAILR1, TNFRSF10A, Tumor necrosis factor receptor superfamily member 10A