



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD261/TNFRSF10A/DR4 Antibody (SAA0083), APC
PRODUCT REFERENCE	ARO-A10497
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	TNF-related apoptosis-inducing ligand receptor 1, APO2, Death receptor 4, TRAIL receptor 1, TRAIL-R1, CD261, DR4, TRAILR1, TNFRSF10A, Tumor necrosis factor receptor superfamily member 10A

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 0.2% BSA, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-TNF-related apoptosis-inducing ligand receptor 1, APO2, Death receptor 4, TRAIL receptor 1, TRAIL-R1, CD261, DR4, TRAILR1, TNFRSF10A, Tumor necrosis factor receptor superfamily member 10A

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY