



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD221/IGF1R Antibody (2.13.2#), FITC

OVERVIEW

Reference	ARO-A10506
Clonality	Monoclonal Antibody
Isotype	IgG2, kappa

PRODUCT INFORMATION

Host Species	Human
Target	CD221, Insulin-like growth factor 1 receptor, IGF1R, Insulin-like growth factor I receptor, IGF-I receptor
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-CD221, Insulin-like growth factor 1 receptor, IGF1R, Insulin-like growth factor I receptor, IGF-I receptor