



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD74 Antibody (SP7219), PE

#### OVERVIEW

Reference	ARO-A10510
Clonality	Monoclonal Antibody
Isotype	IgG1

#### PRODUCT INFORMATION

Host Species	Human
Target	HLA-DR antigens-associated invariant chain, HLA class II histocompatibility antigen gamma chain, CD74, CLIP, Ia antigen-associated invariant chain, DHLAG, Ii
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-HLA-DR antigens-associated invariant chain, HLA class II histocompatibility antigen gamma chain, CD74, CLIP, Ia antigen-associated invariant chain, DHLAG, Ii