



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD18/ITGB2 Antibody (SAA0815), PE

#### OVERVIEW

Reference	ARO-A10511
Clonality	Monoclonal Antibody
Isotype	IgG1

#### PRODUCT INFORMATION

Host Species	Human
Target	Complement receptor C3 subunit beta, Cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta, Integrin beta-2, MFI7, CD18, ITGB2, RHC13101
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Complement receptor C3 subunit beta, Cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta, Integrin beta-2, MFI7, CD18, ITGB2, RHC13101