



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD45RO Rabbit Antibody (UCHL1), PE

OVERVIEW

Reference	ARO-A10537
Clonality	Monoclonal Antibody
Isotype	IgG

PRODUCT INFORMATION

Host Species	Rabbit
Target	T200, Leukocyte common antigen, CD45, L-CA, Receptor-type tyrosine-protein phosphatase C, PTPRC
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-T200, Leukocyte common antigen, CD45, L-CA, Receptor-type tyrosine-protein phosphatase C, PTPRC