



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human CD3E Antibody (UCHT-1), PE

#### OVERVIEW

Reference	ARO-A10572
Clonality	Monoclonal Antibody
Isotype	IgG, kappa

#### PRODUCT INFORMATION

Host Species	Rabbit
Target	T3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T-cell surface glycoprotein CD3 epsilon chain
Form	UCHT1 Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-T3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T-cell surface glycoprotein CD3 epsilon chain