



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Anti-Human IL4 Antibody (SAA0371), PE

#### OVERVIEW

Reference	ARO-A10609
Clonality	Monoclonal Antibody
Isotype	IgG1, kappa

#### PRODUCT INFORMATION

Host Species	Mouse
Target	Binetrakin, B-cell stimulatory factor 1, Pitrakinra, Lymphocyte stimulatory factor 1, IL-4, BSF-1, IL4, Interleukin-4
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-Binetrakin, B-cell stimulatory factor 1, Pitrakinra, Lymphocyte stimulatory factor 1, IL-4, BSF-1, IL4, Interleukin-4