



## PRODUCT DATASHEET

FOR RESEARCH USE ONLY

### Mouse IgG2b, kappa Isotype Control Antibody (MPC-11), PE

#### OVERVIEW

---

Reference	ARO-A10655
Clonality	Monoclonal Antibody
Isotype	IgG2b, kappa

#### PRODUCT INFORMATION

---

Host Species	Mouse
Target	
Form	Liquid
Buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Stability & Storage	4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
Aliases /Synonyms	Anti-0