



PRODUCT DATASHEET

FOR RESEARCH USE ONLY

Anti-Human CD282/TLR2 Antibody (OPN-305), PE

OVERVIEW

| | |
|-----------|---------------------|
| Reference | ARO-A10731 |
| Clonality | Monoclonal Antibody |
| Isotype | IgG4, kappa |

PRODUCT INFORMATION

| | |
|---------------------|---|
| Host Species | Human |
| Target | TIL4, Toll-like receptor 2, TLR2, Toll/interleukin-1 receptor-like protein 4, CD282 |
| Form | Liquid |
| Buffer | 0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300. |
| Stability & Storage | 4°C for short term (1 week) / -20°C or -80°C for long term / an addition of 20-40% glycerol can improves cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| Aliases /Synonyms | Anti-TIL4, Toll-like receptor 2, TLR2, Toll/interleukin-1 receptor-like protein 4, CD282 |