

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human CD254/RANKL/TNFSF11 Antibody (SAA0080), PE
<b>PRODUCT REFERENCE</b>	ARO-A10744
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2, kappa
<b>TARGET</b>	ODF, TNF-related activation-induced cytokine, Osteoclast differentiation factor, TNFSF11, TRANCE, Osteoprotegerin ligand, RANKL, Tumor necrosis factor ligand superfamily member 11, Receptor activator of nuclear factor kappa-B ligand, CD254, OPGL

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Mus musculus (Mouse)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 0.2% BSA, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-ODF, TNF-related activation-induced cytokine, Osteoclast differentiation factor, TNFSF11, TRANCE, Osteoprotegerin ligand, RANKL, Tumor necrosis factor ligand superfamily member 11, Receptor activator of nuclear factor kappa-B ligand, CD254, OPGL